

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27

DoD AFFF01:

DETERMINATION OF PERFLUOROOCCTANOIC ACID
AND PERFLUOROOCCTANESULFONIC ACID IN
AQUEOUS FILM FORMING FOAM (AFFF) FOR
DEMONSTRATION OF COMPLIANCE TO MIL-PRF-24385

Revision 1.0

28 **Authors**

29 Janice L. Willey

30 NAVSEA 04X6 Laboratory Quality and Accreditation Office (LQAO)

31 Janice.willey@navy.mil

32

33 **Acknowledgements**

34 This method was developed under the direction of the Technical Warrant Holder for Fire Protection
35 Systems for the Naval Sea Systems Command (NAVSEA) who expresses appreciation to Dr. Coreen
36 Hamilton and Dr. Bharat Chandramouli of SGS AXYS Analytical Services Ltd for their contributions to
37 this method.

38 The following organizations completed a validation study in their laboratories using this method and
39 provided valuable feedback on the method procedures:

40 Agriculture & Priority Pollutants Laboratories, Inc. (APPL, Inc.) (Clovis, CA)

41 Battelle (Norwell, MA)

42 Bureau Veritas Canada Inc. (Mississauga, Ontario, Canada)

43 General Engineering Laboratories, LLC (GEL, LLC) (Charleston, SC)

44

45 **Disclaimer**

46 Mention of trade names or commercial products does not constitute endorsement or recommendation for
47 use.

48

49

50

51

52

53

54

55

56

57

58 Table of Contents

59 1 Scope and Application..... 4

60 2 Summary of Method..... 5

61 3 Definitions..... 6

62 4 Contamination and Interferences..... 6

63 5 Safety..... 7

64 6 Equipment and Supplies..... 7

65 7 Reagents and Standards..... 9

66 8 Sample Collection, Preservation, Storage, and Holding Times..... 13

67 9 Quality Control..... 16

68 10 Calibration..... 19

69 11 Sample Preparation..... 22

70 12 Instrumental Analysis..... 23

71 13 Data Analysis, Calculation, and Reporting..... 25

72 14 Method Performance..... 27

73 15 Pollution Prevention..... 27

74 16 Waste Management..... 27

75 17 References..... 28

76

77

78

79

80

81

82

83

84

85

86

87

88 1.0 Scope and Application

89 This method describes a procedure for the quantitative determination of two method analytes,
 90 Perfluorooctanoic acid (PFOA) and Perfluorooctanesulfonic acid (PFOS) (Table 1.1), in Class B aqueous
 91 film forming foam (AFFF) concentrates (herein cited as “AFFF samples”). The military performance
 92 specification for Fire Extinguishing Agent, Aqueous Film-Forming Foam (AFFF) Liquid Concentrate, for
 93 Fresh and Sea Water, MIL-PRF-24385 (Reference 17.2), requires the PFOA and PFOS content of AFFF
 94 concentrates be determined (herein cited as “compliance testing”) by a laboratory accredited for this
 95 method under their Department of Defense Environmental Laboratory Accreditation Program (DoD
 96 ELAP) scope of accreditation. Therefore in addition to the requirements contained in this method, the
 97 requirements contained in the latest version of the DoD/DOE Quality Systems Manual for Environmental
 98 Laboratories (DoD/DOE QSM) (Reference 17.1) also apply, with the exception of those contained in
 99 Appendix B, Table B-15. Requirements contained in this method supersede the requirements contained
 100 in Table B-15 for compliance testing purposes.

101
 102 This method utilizes solid phase extraction (SPE) and carbon cleanup techniques to prepare AFFF
 103 samples for analysis by liquid chromatography-tandem mass spectrometry (LC-MS/MS). Selectivity is
 104 optimized through the use of Multiple Reaction Monitoring (MRM) mode and the monitoring of at least
 105 two MS/MS transitions for each method analyte. Isotopically labeled standards of PFOA and PFOS are
 106 used to calibrate and quantify PFOA and PFOS by isotope dilution quantitation. Branched and linear
 107 isomers are included in the quantitation of both method analytes. Inter- and intra-laboratory studies have
 108 generated accuracy and precision data for the determination of PFOA and PFOS in a variety of AFFF
 109 concentrates.

110
 111 This method is restricted to use by or under the supervision of analysts experienced in the use of a liquid
 112 chromatography/tandem quadrupole mass spectrometer and in the interpretation of mass spectra. In
 113 addition, PFAS analysis requires specific skills and experience in minimizing laboratory background and
 114 contamination and it is highly recommended that this method only be performed in laboratories with prior
 115 experience in PFAS analysis by LC-MS/MS.

116 Table 1.1 Method Analyte

Table 1.1 Names, Abbreviations, and CAS Registry Numbers for Method Analytes		
Method Analyte Name	Abbreviation	CAS Number
Perfluorooctanoic acid	PFOA	335-67-1
Perfluorooctanesulfonic acid	PFOS	1763-23-1

117

118 1.1 Detection of PFAS Isomers

119 Per- and polyfluoroalkyl substances (PFAS), including PFOA and PFOS, may exist as branched as well
 120 as linear isomers in AFFF samples, therefore both branched and linear isomers must be included in the
 121 determination of each method analyte. A quantitative standard that contains a mixture of branched and
 122 linear isomers of PFOS is used for calibration and quantification of PFOS. This standard must be used
 123 for all calibration, calibration verifications, and quality control (QC) samples. No such quantitative
 124 standard is currently commercially available for PFOA, therefore a quantitative standard that contains
 125 only the linear isomer of PFOA is used for calibration. A qualitative standard containing a mixture of
 126 branched and linear isomers of PFOA is analyzed post calibration in order to determine the retention time
 127 of the branched isomers of PFOA, to therefore be included in the quantification of PFOA.

128 **1.2 Limit of Quantitation**

129 Single-laboratory and interlaboratory validation studies of this method demonstrated a limit of
 130 quantitation (LOQ) of less than 25 ppb for both PFOA and PFOS is routinely achievable. An LOQ of
 131 less than 25 ppb must be achieved for both PFOA and PFOS in order for results to be applicable to
 132 compliance testing.

133
 134 **1.3 Method Flexibility**

135 The type of LC system (UPLC, HPLC), the LC columns, LC conditions, and MS conditions utilized may
 136 be different than those utilized in the development of this method. Changes may not be made to the
 137 quality control (QC) (Section 9.0), Calibration (Section 10.0), Sample Preparation (Section 11.0),
 138 Instrumental Analysis (Section 12.0), and the Data Analysis, Calculations, and Reporting (Section 13.0)
 139 requirements.

140

141 **2.0 Summary of Method**

142 A 0.02 mL aliquot of the AFFF sample is diluted using PFAS-free reagent water, spiked with an extracted
 143 internal standard (EIS) solution containing isotopically labeled PFOA and PFOS compounds (Table 2.1),
 144 extracted using SPE, and matrix interferences are reduced using a carbon clean-up procedure. The extract
 145 is spiked with non-extracted internal standard (NIS) solution containing other isotopically labeled PFOA
 146 and PFOS compounds (Table 2.1) and analysis is conducted by liquid chromatography-tandem mass
 147 spectrometry (LC-MS/MS) utilizing negative ion spray and multiple reaction monitoring (MRM) mode.

148 This method requires PFOA and PFOS to be quantified and reported in their acid form. The quantitation
 149 scheme utilized by this method is isotope dilution quantitation, which recovery corrects the results for
 150 method analytes using the response of the isotopically labeled PFOA and PFOS compounds (EIS) added
 151 to the sample prior to extraction. To assess overall analytical quality, the recovery of the EIS compound
 152 is determined through comparison of its response to response of the applicable NIS compound that was
 153 added post extraction (see Table 2.2), prior to sample analysis. At a minimum, the transitions listed in
 154 Table 2.2 must be utilized for quantitation and confirmation. A third transition may be added for
 155 confirmation purposes.

156

EIS Compounds	Abbreviation
Perfluoro-n-[¹³ C ₈]octanoic acid	¹³ C ₈ -PFOA
Perfluoro-1-[¹³ C ₈]octanesulfonic acid	¹³ C ₈ -PFOS
NIS Compounds	Abbreviation
Perfluoro-n-[1,2,3,4- ¹³ C ₄]octanoic acid	¹³ C ₄ -PFOA
Perfluoro-n-[1,2,3,4- ¹³ C ₄]octanesulfonic acid	¹³ C ₄ -PFOS

Abbreviation	Example Retention Time	Parent Ion Mass	Quantification Ion Mass	Confirmation Ion Mass	Typical Ion Ratio	Quantification Reference Compound
Analytes						
PFOA	6.16	413.0	369.0	169.0	3.0	¹³ C ₈ -PFOA
PFOS	7.59	498.9	79.9	98.8	2.3	¹³ C ₈ -PFOS
Extracted Internal Standards						

Table 2.2. Analyte Ions Monitored, Extracted Internal Standard, and Non-extracted Internal Standard Used for Quantification						
Abbreviation	Example Retention Time	Parent Ion Mass	Quantification Ion Mass	Confirmation Ion Mass	Typical Ion Ratio	Quantification Reference Compound
¹³ C ₈ -PFOA	6.16	421.1	376.0	NA		¹³ C ₄ -PFOA
¹³ C ₈ -PFOS	7.59	507.1	98.9	79.9		¹³ C ₄ -PFOS
Non-Extracted Internal Standards						
¹³ C ₄ -PFOA	6.16	417.1	172.0	NA		
¹³ C ₄ -PFOS	7.59	502.8	79.9	98.9		

157

158 3.0 Definitions

159 Definitions of terms used in this method are consistent with those published in the DoD/DOE QSM
160 (current version). Additions and deviations from these are provided below.

161 3.1 AFFF

162 Aqueous film-forming foam containing fluorinated and hydrocarbon surfactants and other constituents
163 formulated for class B hydrocarbon fires.

164 3.2 Analysis Batch

165 A set of samples that are analyzed on the same instrument during a 24-hour period.

166 3.3 AFFF Concentrate

167 Formulation of AFFF in concentrate form prior to dilution for end-use firefighting application.

168 3.4 Military Specification

169 A United States defense standard used to help achieve standardization objectives by the U.S. Department
170 of Defense.

171 3.5 PFAS-free

172 Concentration of PFOA and PFOS each have been verified to be less than ½ the Limit of Quantitation
173 (LOQ). This definition is applicable to materials, supplies, and quality control samples (e.g., instrument
174 and method blanks) contained in this method.

175 3.6 Qualitative Standards

176 A standard that contains a mixture of branched and linear isomers of a target analyte that is not of the
177 purity needed to be considered quantitative. These standards are often identified as technical grade
178 mixtures by manufacturers. Qualitative standards are analyzed for comparison to suspected branched
179 isomer peaks in test samples.

180

181

182 4.0 Contamination and Interferences

183 4.1 Labware, Reagents, Equipment, and Supplies

184 Materials, and reagents used and components in the analytical instruments should be handled and chosen
185 carefully to prevent contamination with the method analytes and interferences. Documented best
186 practices to accomplish this are available. Samples in aqueous-based solvents should minimize contact
187 with polypropylene to minimize any adsorption effects. Sample containers must be high density
188 polyethylene (HDPE) containers. The use of glassware for transfer procedures and standards preparation
189 is allowed, but standards and sample extracts are to be stored in HDPE containers.

190 **4.1.1** All glass equipment that is used in the preparation or storage of reagents must be appropriately
191 cleaned prior to use. The process must include detergent, reagent water, and solvent rinses
192 followed by baking in a kiln or furnace. After detergent washing, glassware should be rinsed
193 immediately with reagent water. A solvent rinse procedure using methanolic ammonium
194 hydroxide (1%), toluene, and methanol is recommended.

195 **4.1.2** Due to the potential for sample preparation and analysis supplies to introduce PFAS into the
196 sample, the residual PFAS content of disposable plasticware and filters must be verified by
197 batch/lot number. If the residual PFAS content of these materials is at a concentration less than
198 half the limit of quantitation, these materials can be utilized without cleaning.

199 **4.1.3** The SPE manifold can be a significant source of PFAS contamination. All parts of the SPE
200 manifold must be cleaned between samples by sonicating in methanolic ammonium hydroxide
201 (1%) and air drying prior to use. In order to validate the cleaning process used is sufficient, the
202 manifold position of the batch quality control (QC) samples (specifically the method blank (MB)
203 and laboratory control sample (LCS) in each batch must be rotated on a batch by batch basis.

204 **4.2** Samples

205 **4.2.1** AFFF samples can contain high concentrations of PFAS, therefore the risk of PFAS cross-
206 contamination is high. To reduce the potential for cross-contamination, all equipment must be
207 cleaned prior to, and after each use. The cleaning solvents typically used include water,
208 methanol, and methanolic ammonium hydroxide.

209 **4.2.2** Interferences co-extracted from samples will vary considerably sample to sample, depending on
210 the non-PFAS constituents of each AFFF concentrate. Concentrations of these interfering
211 compounds can be several orders of magnitude higher than concentration of the targeted PFAS
212 analytes. Given the goal of this method is to achieve quantitation of low levels of PFAS, it is
213 critical that these interferences are eliminated to the greatest extent possible. The SPE and carbon
214 cleanup procedures contained in this method is included for this purpose.

215 **5.0** Safety

217 **5.1** The toxicity or carcinogenicity of all of the PFAS included in the scope of this method are yet to
218 be determined; therefore they should be treated as a potential health hazard. Requirements contained in
219 the laboratory's Safety Manual and associated standard operating procedures (SOPs) should be followed.

220 **5.2** A safety data sheets (SDS) must accompany each AFFF sample submitted and be retained by the
221 laboratory in accordance with Occupational Safety and Health Administration (OSHA) regulations and
222 laboratory policy.

223

224 **6.0** Equipment and Supplies

- 225 Brand names, suppliers, and part numbers are for illustration purposes only and no endorsement is
226 implied. If equivalent performance can be achieved using apparatus and materials other than those
227 specified here, they can be used.
- 228 **6.1** Sample bottles and caps –HDPE, with liner-less HDPE or polypropylene caps. Use of PTFE-
229 lined caps is prohibited.
- 230 **6.2** Nitrile gloves
- 231 **6.3** Balances
- 232 **6.3.1** Analytical – Capable of weighing 0.1 mg
- 233 **6.3.2** Top loading – Capable of weighing 10 mg
- 234 **6.4** Ultrasonic mixer (sonicator)
- 235 **6.5** HDPE bottles, with liner-less HDPE or polypropylene caps – 60 mL and 500 mL
- 236 **6.6** pH Paper, range 0-14 - (Whatman® Panpeha™ or equivalent), 0.5-unit readability
- 237 **6.7** Analog or digital vortex mixer, single or multi-tube (Fisher Scientific 02-215-452, or equivalent)
- 238 **6.8** Volumetric flasks, Class A
- 239 **6.9** Disposable polypropylene collection tubes (13 x 100 mm, 8 mL)
- 240 **6.10** Variable speed mixing table (Fisherbrand™ Nutating mixer or equivalent)
- 241 **6.11** Silanized glass wool (Sigma-Aldrich, Cat # 20411 or equivalent) – store in a clean glass jar and
242 rinsed with methanol (2 times) prior to use.
- 243 **6.12** Disposable syringe filter, 25-mm, 0.2-µm Nylon membrane, PALL/Acrodisc or equivalent
- 244 **6.13** Glass fiber filter, 47 mm, 1 µm, PALL A/E or equivalent
- 245 **6.14** Centrifuge (Thermo Scientific Legend RT+, 16 cm rotor, or equivalent), capable of reaching at
246 least 3000 rpm
- 247 **6.15** Centrifuge tubes – Disposable HDPE or polypropylene centrifuge tubes (15 and 500 mL)
- 248 **6.16** Norm-Ject® syringe (or equivalent), polypropylene/HDPE, 5 mL
- 249 **6.17** Variable volume pipettes with disposable HDPE or polypropylene tips (10 µL to 5 mL) – used for
250 preparation of calibration standards and spiked samples.
- 251 **6.18** Disposable glass pipets
- 252 **6.19** Calibrated mechanical pipettes or Hamilton graduated syringes
- 253 **6.20** Solid-phase extraction (SPE) cartridges (Waters Oasis WAX 150 mg, Cat # 186002493 or
254 equivalent). The SPE sorbent must have a pKa above 8 so that it remains positively charged
255 during the extraction.
- 256 **6.21** Vacuum manifold for SPE Cartridges (Waters™ extraction manifold #WAT200607 or
257 equivalent)

- 258 **6.22** Automatic or manual solvent evaporation system (TurboVap® LV or equivalent)
- 259 **6.23** Evaporation/concentrator tubes: 60 mL clear glass vial, 30 x 125 mm, without caps
260 (Wheaton Cat # W226060 or equivalent). Cover with foil if required.
- 261 **6.24** Snap cap/crimp top vials, 300 µL, polypropylene (12 x 32 mm) – used in sample pre-
262 screening (DWK Life Sciences Cat # 225180 or equivalent)
- 263 **6.25** Polypropylene crimp/snap vials, 1 mL (Agilent Cat # 5182-0567 or equivalent)
- 264 **6.26** Clear snap cap, PVDC film/white silicone, 11 mm (American Chromatography
265 Supplies Cat # C299-11 or equivalent)
- 266 **6.27** Single step filter vials (Restek Thomson SINGLE StEP® Standard Filter Vials, 0.2 µm
267 Nylon membrane, with Black Preslit caps Cat # 25891 or equivalent) – used in sample
268 pre-screening
- 269 **6.28** 10 mg polypropylene or stainless steel scoops
- 270 **6.29** Ultra high-performance liquid chromatograph (UPLC also called UHPLC) or high- performance
271 liquid chromatograph (HPLC) equipped with tandem quadrupole mass spectrometer (Waters
272 Xevo TQ-S Micro or equivalent capable of electrospray ionization in the negative ion mode.
- 273 **6.30** C18 column, 1.7 µm, 50 x 2.1 mm (Waters Acquity UPLC® BEH or equivalent)
- 274
- 275 **6.31** Guard column (Phenomenex Kinetex® Evo C18 or equivalent)
- 276
- 277 **6.32** Trap/delay column (Purospher Star RP-18 endcapped [3 µm] Hibar® RT 50-4 or equivalent)
- 278
- 279 **6.33** Bottles, HDPE or glass, with liner-less HDPE or polypropylene caps. Various sizes. To store
280 prepared reagents.
- 281
- 282 **7.0 Reagents and Standards**
- 283 **7.1 Reagents**
- 284 When prepared by the laboratory, must be stored in either glass or HDPE containers.
- 285 **7.1.1** Acetic acid - ACS grade or equivalent, store at room temperature
- 286 **7.1.2** Acetic acid (0.1%) - dissolve acetic acid (1 mL) in reagent water (1 L), store at room temperature,
287 replace after 3 months. This reagent is used for sample extract dilution only.
- 288 **7.1.3** Acetonitrile – UPLC grade or equivalent, verified before use, store at room temperature
- 289 **7.1.4** Ammonium acetate - (Caledon Ultra LC/MS grade, or equivalent), store at 2-8° C, replace 2
290 years after opening date
- 291 **7.1.5** Ammonium hydroxide - certified ACS+ grade or equivalent, 30% in water, store at room
292 temperature
- 293 **7.1.6** Aqueous ammonium hydroxide (3%) - add ammonium hydroxide (10 mL, 30%) to reagent water
294 (90 mL), store at room temperature, replace after 3 months

- 295 **7.1.7** Methanolic ammonium hydroxide (0.3%) - add ammonium hydroxide (1 mL, 30%) to methanol
296 (99 mL), store at room temperature, replace after 1 month
- 297 **7.1.8** Methanolic ammonium hydroxide (1%) - add ammonium hydroxide (3.3 mL, 30%) to methanol
298 (96.7 mL), store at room temperature, replace after 1 month
- 299 **7.1.9** Methanolic ammonium hydroxide (2%) - add ammonium hydroxide (6.6 mL, 30%) to methanol
300 (93.4 mL), store at room temperature, replace after 1 month
- 301 **7.1.10** Methanolic potassium hydroxide (0.05 M) – add 3.3 g of potassium hydroxide to 1 L of
302 methanol, store at room temperature, replace after 3 months
- 303 **7.1.11** Methanol with 4% water, 1% ammonium hydroxide and 0.625% acetic acid - add ammonium
304 hydroxide (3.3 mL, 30%), reagent water (1.7 mL) and acetic acid (0.625 mL) to methanol (92
305 mL), store at room temperature, replace after 1 month. This solution is used to prepare the
306 instrument blank (Section 7.3.6).
- 307 **7.1.12** Eluent A – Acetonitrile, Caledon Ultra LCMS grade or equivalent
- 308 **7.1.13** Eluent B - 2 mM ammonium acetate in 95:5 water/acetonitrile. Dissolve 0.154 g of ammonium
309 acetate (Section 7.1.4) in 950 mL of water and 50 mL of acetonitrile (Caledon Ultra LCMS grade,
310 or equivalent). Store at room temperature, shelf life 2 months.
- 311 **7.1.14** Formic acid - (greater than 96% purity or equivalent), verified by lot number before use, store at
312 room temperature
- 313 **7.1.15** Formic acid (aqueous, 0.1 M) - dissolve formic acid (4.6 g) in reagent water (1 L), store at room
314 temperature, replace after 2 years
- 315 **7.1.16** Formic acid (aqueous, 0.3 M) - dissolve formic acid (13.8 g) in reagent water (1 L), store at room
316 temperature, replace after 2 years
- 317 **7.1.17** Formic acid (aqueous, 5% v/v) - mix 5 mL formic acid with 95 mL reagent water, store at room
318 temperature, replace after 2 years
319
- 320 **7.1.18** Formic acid (aqueous, 50% v/v) - mix 50 mL formic acid with 50 mL reagent water, store at
321 room temperature, replace after 2 years
322
- 323 **7.1.19** Formic acid (methanolic 1:1, 0.1 M formic acid/methanol) - mix equal volumes of methanol and
324 0.1 M formic acid, store at room temperature, replace after 2 years
325
- 326 **7.1.20** Methanol - (HPLC grade or better, 99.9% purity), verified by lot number before use, store at
327 room temperature
328
- 329 **7.1.21** Potassium hydroxide – certified ACS or equivalent, store at room temperature, replace after 2
330 years
331
- 332 **7.1.22** Reagent water – Laboratory reagent water, test by lot/batch number for residual PFAS content
333
- 334 **7.1.23** Carbon – EnviCarb® 1-M-USP or equivalent, verified by lot number before use, store at room
335 temperature. Loose carbon allows for better adsorption of interferent organics. The use of
336 carbon cartridges is not allowed.
337

338 **7.1.24** Reference matrix – PFAS-free reagent water, purified water, Type I. Used to prepare the batch
 339 QC samples (e.g., method blank, limit of quantitation verification sample, and laboratory control
 340 samples)

341
 342 **7.2** Standards

343 Solutions are prepared by the laboratory using solutions or mixtures (prime stocks) with certification to
 344 their purity, concentration, and authenticity. Standard solutions must be stored in the dark at less than 4
 345 °C unless the vendor recommends otherwise in screw-capped vials with foiled-lined caps. Monitor
 346 solutions for evaporation. If loss is detected, replace the solution.

347 **7.2.1** Native Spiking Standard Solutions

348 Prepare the native spiking standard by diluting prime stocks that contain PFOA and PFOS with
 349 methanol. Quantitative prime stock solutions must be used to create this standard. As stated in
 350 Section 1.1, the quantitative standard used for PFOS must contain a mixture of branched and
 351 linear isomers while the quantitative standard used for PFOA contains only the linear isomer. It is
 352 used to prepare the calibration, instrument sensitivity check, initial calibration verification and
 353 continuing calibration verification standards and to spike preparation batch QC samples (LCS,
 354 LCSD, and LOQVER).

355 **7.2.2** Qualitative Standard

356 As stated in Section 1.1, a qualitative standard that contains a mixture of branched and linear
 357 isomers of PFOA must be analyzed prior to sample analysis. Prepare this standard by diluting a
 358 qualitative stock standard containing PFOA with a solution that matches the solvent mix of
 359 sample extracts, which contain methanol with 4% water, 1% ammonium hydroxide and 0.625%
 360 acetic acid. This standard is used for comparison with suspected branched isomer peaks in AFFF
 361 samples.

362 **7.2.3** Extracted Internal Standards (EIS) Solutions

363 Prepare extracted internal standard solutions by diluting, with methanol, prime stock standards
 364 containing the isotopically labeled compounds listed in Table 7.1. Table 7.1 provides the volume
 365 of EIS solution used to spike samples and the resulting nominal amount of each compound in the
 366 sample.

EIS Compounds	Volume Spiked	Amount Added (ng)
¹³ C ₈ -PFOA	50 µL	10
¹³ C ₈ -PFOS	50 µL	10

367
 368 **7.2.4** Non-Extracted Internal Standard (NIS) Solutions

369 Prepare non-extracted internal standard solutions by diluting, with methanol, prime stock
 370 standards containing the isotopically labeled compounds listed in Table 7.2. Table 7.2 provides
 371 the volume of NIS solution used to spike samples and the resulting nominal amount of each
 372 compound in the sample.

NIS Compounds	Volume Spiked	Amount Added (ng)
¹³ C ₄ -PFOA	50 µL	10
¹³ C ₄ -PFOS	50 µL	10

373

374 **7.2.5 Calibration Standard Solutions**

375 Prepare a minimum of 6 calibration standard solutions by diluting native standards with
 376 methanol, methanolic ammonium hydroxide (2%), water, and acetic acid to achieve final
 377 concentrations of PFOA and PFOS that encompass the working range of the instrument. NIS and
 378 EIS compounds are added to each calibration standard such that the concentration of these
 379 compounds remain constant over the series of calibration standards. After dilution, the final
 380 solutions will match the solvent mix of sample extracts, which contain methanol with 4% water,
 381 1% ammonium hydroxide and 0.625% acetic acid. Calibration standard solutions do not undergo
 382 solid phase extraction or cleanup.

383 A minimum of five contiguous calibrations standards are required for a valid analysis when using
 384 a linear calibration model, with at least five calibration standards within the quantitation range
 385 (i.e., from the LOQ to the highest calibration standard). If a second-order calibration model is
 386 used, then a minimum of six calibration standards are required, with at least six calibration
 387 standards within the quantitation range. The lowest level calibration standard must meet a signal-
 388 to-noise ratio of 3:1 and be at a concentration less than or equal to the Limit of Quantitation
 389 (LOQ). Table 7.3 provides the concentrations for the eight calibration solutions utilized during
 390 method development. If instrument sensitivity allows, additional lower concentration standards
 391 may be added to accommodate a lower LOQ.
 392

Analytes	CS1	CS2	CS3	CS4	CS5(CC V)	CS6	CS7	CS8
PFOA	0.05	0.1	0.2	0.5	2.5	12.5	62.5	250
PFOS	0.05	0.1	0.2	0.5	2.5	12.5	62.5	250
Extracted Internal Standard (EIS) Compounds								
¹³ C ₈ -PFOA	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
¹³ C ₈ -PFOS	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Non-extracted Internal Standard (NIS) Compounds								
¹³ C ₄ -PFOA	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
¹³ C ₄ -PFOS	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5

393

394 **7.2.6 Instrument Sensitivity Check (ISC) Standard**

395 Prepare an instrument sensitivity check (ISC) standard in the same manner as the calibration
396 standards solutions. The concentration of PFOA and PFOS in the ISC standard must correspond
397 to the LOQ for each method analyte.

398 **7.2.7** Initial Calibration Verification (ICV) Standard

399 Initial calibration verification standards are to be prepared in the same manner as calibration
400 standards, however, must use a prime stock for a source different than the prime stock used to
401 prepare the calibration standards. This alternative source can be either a different manufacturer
402 or a different lot from the same manufacturer. The concentration of PFOA and PFOS can range
403 from the LOQ to the mid-level calibration standard concentration.

404 **7.2.8** Continuing Calibration Verification (CCV) Standard

405 Prepare the continuing calibration verification standards in the same manner as the calibration
406 standards. The concentration of PFOA and PFOS can range from the LOQ to the mid-level
407 calibration standard concentration. The mid-level calibration standard solution can serve as the
408 CCV.

409 **7.2.9** Instrument Blank

410 Prepare instrument blanks (IBs) by fortifying the solution in Section 7.1.11 with EIS and NIS.

411

412 **8.0 Sample Collection, Preservation, Storage, and Holding Times**

413 **8.1** Scope of Sampling Procedure

414 This sample collection procedure is applicable to sampling performed by the AFFF concentrate vendor
415 prior to submission of the concentrate for consideration, by Naval Research Laboratory (NRL) personnel
416 during qualification inspection, conformance inspection, retention inspection, and by facilities for
417 verification of composition of the formulation of the lot purchased is valid only when it conforms to the
418 requirements herein. All samples, regardless of the entity collecting the sample, must be collected from
419 containers that are compliant with MIL-PRF-24385, section 3.6 requirements.

420 It is critical that the sampling of current and future AFFF stock concentrates be performed in a safe and
421 consistent manner such that the subsequent analysis is representative of the original concentrate and the
422 concentrations of the constituents within each product meet product specifications.

423 **8.2** Sampling Procedure Summary

424 Prior to sampling, drums or containers must be inspected and opened by personnel possessing a thorough
425 understanding of MIL-PRF-24385 requirements and PFAS sampling precautions and practices.

426 Inspection involves the observation and recording of visual qualities of each drum/container and any
427 characteristics or identification markings pertinent to the classification of the drum's contents. Full
428 documentation is captured on a chain-of-custody (COC). Sampling should be performed on previously
429 unopened AFFF stock concentrate containers if possible, in order to avoid any question of
430 providence/integrity. Sampling of a small but representative volume of the container is to be performed
431 using a glass drum thief sampler and the contents placed into a high density polyethylene (HDPE) bottle
432 with a polypropylene, unlined cap.

433 8.3 AFFF Concentrate Stock Containers and Sample Bottles

434 Sample bottles (and lids) for collection of concentrated AFFF must be virgin HDPE bottles verified to be
435 PFAS-free and supplied by the laboratory performing the analysis. HDPE bottles with unlined
436 polypropylene caps are the only bottles that may be used for these samples. Use of any other types of
437 bottles, especially Teflon™ sample bottles or Teflon-lined lids of any kind will void the analysis and
438 require resampling. Avoid agitating the AFFF drum prior to sampling in order to avoid foaming. AFFF
439 drums should be upright and relatively level. The drum must be easily identifiable and preferably
440 unopened.

441 8.4 Precautions

442 All personnel should review the SDS of the specific AFFF concentrate being sampled and be warned of
443 the hazards prior to handling AFFF drums. An adequate volume of absorbent material should be kept near
444 areas where minor spills may occur. Where major spills may occur, a containment berm adequate to
445 contain the entire volume of liquid in the drums should be constructed before any handling takes place. If
446 drum contents spill, personnel trained in spill response should be used to isolate and contain the spill.

447 8.5 Equipment and Supplies

448 **8.5.1** Personal protection equipment (nitrile gloves, goggles, etc.)

449 **8.5.2** Virgin 75-150 mL, PFAS-free HDPE sample bottle with un-lined screw cap

450 **8.5.3** Chain of Custody

451 **8.5.4** Virgin glass drum thief sampler (no larger than 12 mm O.D.)

452 **8.5.5** Drum opening devices

453 **8.5.6** Sample Labels

454 **8.5.7** Sampling notebook

455 **8.5.8** Re-sealable Plastic Bags

456 **8.5.9** Paper towels

457 **8.5.10** Sample bottle custody seals

458 8.6 Preparation and Inspection

459 AFFF containers/drums should be visually inspected to gain as much information as possible about their
460 contents. The drums should be inspected for general condition, punctures, leaking contents, evidence of
461 broken seals and signs that the drum is under pressure. It is recommended, but not required, that
462 photographs of the drum, including the drum labels and drum opening (prior to breaking the seal) be
463 taken and kept as part of the sampling record. Ensure the equipment and supplies listed in Section 8.5 are
464 readily available and quantities are sufficient for the number of drums to be sampled.

465 8.7 Sampling Procedure

466 Since some layering or stratification may occur in any solution left undisturbed and mixing could cause
467 excessive foaming, a drum thief is used to obtain a sample that represents the entire depth of the
468 container. A glass thief sampler is most widely used for sampling container liquids. It is a simple glass
469 tube and is commonly referred to as a drum thief. This tool is cost effective, quick, and disposable. A new
470 drum thief **MUST** be used for each container in order to decrease the chance for cross contamination. A
471 drum thief that is used for AFFF sampling is not reusable and must be disposed of after sampling.
472 Appropriate personal protective equipment such as nitrile gloves and goggles, at a minimum, must be
473 worn during the sampling process. Follow steps 8.7.1 through 8.7.13 for each AFFF container to be
474 sampled.

475 **8.7.1** Record all identifying information from the AFFF container label including the product name,
476 manufacturer, batch number, date of manufacture, location of manufacture, expiration date
477 (month/year), if applicable, concentrate type (3% or 6%), and the condition of the seal in the
478 sampling notebook.

479 **8.7.2** Create a unique identifier that will be used to track the sample through the sampling and analysis
480 process. Record this identifier, the date of sampling and the sampler's initials in the sampling
481 notebook, on the COC, on a label, and on the original AFFF container. It is recommended, but not
482 required, that photographs be taken of the AFFF container clearly showing all information on the
483 label and the unique identifier and of the opening with the seal intact be taken to further
484 document the sampling event.

485 **8.7.3** Break the seal and open the container very slowly to allow for the gradual release of any built-up
486 pressure.

487 **8.7.4** Remove the cap from one of the laboratory supplied sample bottles.

488 **8.7.5** Slowly insert drum thief as deep into the drum as possible, keeping the drum thief in a completely
489 vertical orientation and allow the AFFF in the container to reach natural level in the drum thief.

490 **8.7.6** Cap the top of the drum thief with a tapered stopper or thumb; ensuring liquid does not come into
491 contact with stopper.

492 **8.7.7** Carefully remove the capped drum thief from the AFFF container. Then insert the bottom of the
493 drum thief into the sample bottle.

494 **8.7.8** Release thumb or stopper and allow the drum thief to completely drain into the sample bottle.
495 Repeat as necessary until the bottle is approximately two-thirds full. Do not allow a subsample to
496 overflow the sample bottle.

497 **8.7.9** Dispose of glass drum thief according to the appropriate procedures specified by the facility
498 management.

499 **8.7.10** Reseal the AFFF drum.

500 **8.7.11** Cap the sample bottle tightly and wipe off any excess AFFF from the outside of the bottle with
501 clean, unused paper towels.

502 **8.7.12** Record the time of sampling in the sampling notebook, on the COC, on the drum, and on the
503 label. Place the label on the sample bottle.

504 **8.7.13** Place a custody seal over the cap/bottle interface and then place the sample bottle into a plastic
505 bag and seal.

506 **8.7.14** After all of the AFFF samples have been taken, ensure all other appropriate information has been
507 recorded on the COC, the information recorded on the COC matches the information on the
508 corresponding sample bottle, and the COC has been signed by the personnel performing the
509 sampling. The sampler must retain a copy of the completed COC; however, the original must be
510 included in the shipping container and follow the sample to the laboratory. Package the samples
511 for transport per the laboratory's instructions.

512 **8.8** Sample Preservation, Shipping, and Holding Times

513 **8.8.1** Sample Preservation and Shipping

514 AFFF samples collected from containers are in concentrated form and do not require any
515 chemical preservation. In addition, these concentrated samples are not subject to excessive
516 thermal degradation or decomposition and once sealed can be kept and transported without
517 thermal preservation. There is currently limited data on whether the PFAS constituents are
518 photosensitive therefore, the sample should be protected from direct sunlight.

519 **8.8.2** Holding Times

520 No formal Holding Time Studies of PFAS content in AFFF concentrates has been published to
521 date. This method takes a conservative approach, requiring collected samples be prepared within
522 90 days of collection. AFFF sample extracts must be stored in the dark at less than 4 °C until
523 analyzed. If stored in the dark at less than 4 °C, sample extracts may be stored for up to 30 days
524 prior to analysis.

525

526 **9.0 Quality Control**

527 **9.1** Initial Demonstration of Capability

528 DoD ELAP accreditation for the analysis of PFOA and PFOS in AFFF concentrates in accordance with
529 this method is required for PFOA and PFOS content compliance testing per MIL-PRF-24385. As such,
530 laboratories performing this method must meet the requirements of the DoD/DOE QSM (current version).
531 These requirements include, but are not limited to, those applying to Initial Demonstration of Capability
532 (IDC), detection limit (DL) determination, Limit of Detection (LOD) verification, LOQ verification,
533 reporting requirements, and proficiency testing. The requirements of this method supersede the PFAS-
534 specific ongoing QC requirements contained in the DoD/DOE QSM, Version 5.3, Appendix B, Table B-
535 15.

536 **9.2** Ongoing QC Requirements

537 This section describes the ongoing QC elements that must be included when processing and analyzing
538 field samples. Table A.1 in Appendix A provides a summary of the acceptance criteria for each QC
539 sample.

540 9.2.1 Analytical Batch QC Samples**541 9.2.1.1 Instrument Blank (IB)**

542 IBs (Section 7.2.9) evaluate the background concentrations of method analytes associated
543 with the analytical system. Background concentrations of method analytes must be less
544 than one-half the LOQ. EIS compounds must recover within 50–200% of their true
545 value, and NIS compounds must recover within 50–200% of their true value. Instrument
546 blanks (IBs) (Section 7.2.9) must be analyzed before analysis of the calibration curve,
547 after the highest-level calibration standard, and after CCVs. In addition, IB(s) must be
548 analyzed following samples whose PFOA and/or PFOS concentrations exceed the
549 quantitation range, until method analyte concentrations are less than or equal to one-half
550 the LOQ in an IB. Any samples analyzed before this criteria is achieved, must be re-
551 analyzed using a new aliquot of the final sample extract, if the exceeding method
552 analyte(s) is at a concentration greater than or equal to the LOQ. Subtraction of IB
553 concentrations from sample results is not permitted.

554 9.2.1.2 Instrument Calibration Verification (ICV)

555 An ICV (Section 7.2.7) must be analyzed after the analysis of the IB which follows the
556 calibration curve. The acceptance criterion for the ICV is 70–130% of the true value of
557 PFOA and PFOS, EIS compounds must recover within 50–200% of their true value, and
558 NIS compounds must recover within 50–200% of their true value. If this criteria are not
559 met, a fresh ICV should be prepared and analyzed. Analysis of samples cannot proceed
560 until the analysis of an ICV meeting this criterion.

561 9.2.1.3 Instrument Sensitivity Check (ISC)

562 An ISC (Section 7.2.6) must be analyzed after each calibration curve and daily, following
563 the first IB of the analytical sequence. The acceptance criteria for the ISC is 70–130% of
564 the true value of PFOA and PFOS EIS compounds must recover within 50–200% of their
565 true value, and NIS compounds must recover within 50–200% of their true value. If these
566 criteria are not met, a fresh ISC should be prepared and analyzed. Analysis of samples
567 cannot proceed until the analysis of an ISC meeting this criterion.

568 9.2.1.4 Continuing Calibration Verification (CCV)

569 A CCV (Section 7.2.8) must be analyzed at the beginning of each analysis batch, after
570 every ten field samples, and at the end of the analysis batch. The acceptance criterion for
571 the CCV is 70–130% of the true value of PFOA and PFOS EIS compounds must recover
572 within 50–200% of their true value, and NIS compounds must recover within 50–200%
573 of their true value. If a CCV exceeds this criterion, immediately analyze two additional
574 consecutive CCVs. If both of these CCVs are within the criteria, samples may be
575 reported without re-analysis. If either exceed the criteria, or if two consecutive CCVs
576 cannot be analyzed immediately after the failing CCV, corrective action must be taken.
577 Once correction has been made and a CCV has been analyzed and has met the criteria, all
578 samples bracketed by the failing CCV must be reanalyzed.

579 9.2.2 Preparation Batch QC Samples

580 A preparation batch consists of up to 20 AFFF samples (10 unique AFFF samples and their
581 associated duplicates) and preparation batch QC samples that are extracted together using the
582 same lot of each material used (e.g., solid phase cartridge, fortifying solutions, and solvents). The
583 preparation batch QC samples included in each batch must consist of, at a minimum, a Method
584 Blank (MB), Laboratory Control Sample (LCS), Laboratory Control Sample Duplicate (LCSD),
585 and LOQ Verification (LOQVER), and duplicate samples for each AFFF sample included in the
586 batch. PFAS-free reagent water is used as the reference media for the MB, LCS, LCSD, and
587 LOQVER.

588 **9.2.2.1 Method Blank (MB)**

589 The concentration of PFOA and PFOS each must be < one-half the LOQ. If the
590 concentration of PFOA and/or PFOS in the MB exceeds this criterion and the
591 concentration of that analyte(s) is equal to or greater than the LOQ in the AFFF sample,
592 the corresponding AFFF result is not valid and the AFFF sample must be re-extracted.
593 Subtracting blank values from sample results is not permitted.

594 **9.2.2.2 Laboratory Control Sample/Laboratory Control Sample Duplicate (LCS/LCSD)**

595 The recoveries of PFOA and PFOS in the LCS and LCSD must be within 70–130% of
596 their true value. The relative percent difference (RPD) for both the PFOA and PFOS
597 results for the LCS and LCSD must be $\leq 30\%$. If these criteria are not met, then all
598 samples in the associated batch must be re-extracted.

599 **9.2.2.3 Limit of Quantitation (LOQ) Verification Sample**

600 The recoveries of PFOA and PFOS in the LCS and LCSD must be within 70–130% of
601 their true value. If these criteria are not met, then all samples in the associated batch
602 must be re-extracted.

603 **9.2.2.4 Sample Duplicate (SD)**

604 For concentration of PFOA and/or PFOS that are equal to or greater than the LOQ, the
605 relative percent difference (RPD) between the sample and the corresponding SD must be
606 $\leq 30\%$. If this criterion is exceeded, then the sample and associated SD must be re-
607 extracted.

608 **9.2.2.5 Extracted Internal Standard (EIS) Recoveries**

609 The concentration of Extracted Internal Standard (EIS) compounds are quantitated with
610 respect to the non-extracted internal standard (NIS) response. EIS compounds must
611 recover within 50–200% of their true value. If the recovery of EIS compounds falls
612 outside of this range, the sample and associated SD must be re-extracted. If the
613 exceedance is associated with the MB, LCS, LCSD, and/or LOQVER the entire batch
614 must be extracted. Samples bracketed by a calibration verification (ICV, CCV, or ISC)
615 that fail to meet these criteria must be reanalyzed.

616 **9.2.2.6 Non-extracted Internal Standard (NIS) Recoveries**

617 The recovery of the Non-extracted Internal Standard (NIS) compounds is determined by
618 comparing the NIS compound peak areas with the average area of the corresponding NIS
619 in the calibration curve standards. NIS compounds must recover within 50–200% of this

620 average. If the recovery of a NIS compound falls outside of this range, the sample (and
 621 associated SD) must be re-extracted. If the exceedance is associated with the MB, LCS,
 622 LCSD, and/or LOQVER the entire batch must be extracted. Samples bracketed by a
 623 calibration verification (ICV, CCV, or ISC) that fail to meet these criteria must be
 624 reanalyzed.

625

626 10.0 Calibration

627 10.1 Mass Calibration and Mass Calibration Verification

628 **10.1.1** Mass calibration must be performed at least annually and must be repeated on an as-needed basis
 629 (e.g., QC failures, ion masses fall outside of the instrument required mass window, major
 630 instrument maintenance, or if the instrument is moved) in accordance with procedures prescribed
 631 by the manufacturer. The manufacturer's instructions for confirmation of the mass calibration,
 632 mass resolution, and peak relative response must be followed and all criteria must be met. The
 633 procedures used must evaluate an ion range that encompasses the ion range (Q1 and Q2 m/z) of
 634 the analytes of interest and isotopically labeled compounds of this method (Table 2.2).

635 **10.1.2** While the MS conditions used during the development of this method (Table 10.1) are provided
 636 below for guidance, suitable operating conditions must be established in accordance with the
 637 manufacturer's instructions.

MS/MS Conditions	Source Temp (°C)	140
	Desolvation Temp (°C)	500
	Capillary Voltage (kV)	0.70
	Cone Gas (L/h)	~70
	Desolvation gas (L/h)	~800

638

639 **10.1.3** Each mass calibration performed must be followed by a mass calibration verification prior to
 640 standards and samples analysis. These verifications must be performed in accordance with
 641 manufacturer's instructions. The ion masses evaluated to the manufacturer's acceptance criteria
 642 must bracket the masses of interest (quantitative and qualitative ions).

643 10.2 Chromatographic Conditions

644 **10.2.1** Chromatographic conditions must be optimized for compound separation and sensitivity.
 645 Conditions must be kept the same for the analysis of all standards, blanks, QC samples, and
 646 AFFF samples. Table 10.2 provides the instrumentation and chromatographic conditions utilized
 647 in the development of this method. Operating conditions are dependent on the instrumentation
 648 used. The LC gradient program runtime of 12.0 minutes helps to ensure later eluting components
 649 found in AFFF samples do not adversely affect the superseding sample analyses. This capability
 650 must be demonstrated for the LC gradient program used as part of the optimization process.

Table 10.2 Instrumentation and General LC Operating Conditions			
Instrument	Waters Acquity UPLC, TQ-S Xevo MS/MS or TQ-S Xevo Micro MS/MS, (or equivalent)		
LC Column	Waters Acquity UPLC @ BEH C18 column, 1.7 µm, 50 x 2.1 mm (or equivalent)		
Guard Column	Guard column: Phenomenex Evo C18 Guard (or equivalent)		
Acquisition	MRM mode, negative ESI, unit resolution		
Injection Volume	2.0 µL (Note: Injection volume will vary with instrument, this volume is provided for the default method only.)		
General LC Conditions			
Column Temp (°C)	40		
Max Pressure (bar)	1100.0		
LC Gradient Program			
Time (min)	Flow mixture ^{1,2}	Flow Rate Program	Gradient Curve
0.0	2% eluent A, 98% eluent B	0.35 mL/min	Initial
0.2	2% eluent A, 98% eluent B	0.35 mL/min	2
4.0	30% eluent A, 70% eluent B	0.40 mL/min	7
7	55% eluent A, 45% eluent B	0.40 mL/min	8
9	75% eluent A, 25% eluent B	0.40 mL/min	8
10	95% eluent A, 5% eluent B	0.40 mL/min	6
10.4	2% eluent A, 98% eluent B	0.40 mL/min	10
11.8	2% eluent A, 98% eluent B	0.40 mL/min	7
12.0	2% eluent A, 98% eluent B	0.35 mL/min	1

¹ Eluent A = Acetonitrile

² Eluent B = 2 mM ammonium acetate in 95:5 water/acetonitrile

651
652

653 10.2.2 Retention Time Calibration

654 **10.2.2.1** The retention time of PFOA, PFOS, and their associated EIS and NIS compounds are
 655 provided in Table 2.2. Since retention time is dependent on the columns/mobile phase
 656 combination employed, the elution order of PFOA, PFOS, and their associated NIS and
 657 EIS compounds must be verified for the combination employed by injecting a series of
 658 solutions, each containing a single compound. On the days an initial calibration (ICAL)
 659 is performed, the retention time of each compound is established by the retention time of
 660 the mid-level standard of the ICAL. On days an ICAL is not performed, the retention
 661 time of the initial CCV of the day or the mid-level standard of the ICAL can be used to
 662 establish the retention time.

663 **10.2.2.2** The retention time of method analytes and NIS and EIS compounds must fall within 0.4
 664 minutes of the established retention time. On days an ICAL is performed, the RT for each
 665 analyte, EIS compound, and NIS compound shall be set using the midpoint standard of
 666 the ICAL curve. On days an ICAL is not performed, the initial daily continuing
 667 calibration verification for that day or the midpoint standard of the ICAL curve can be
 668 used to establish the RT. All branched isomers identified in the PFOS calibration
 669 standards and the qualitative PFOA standard must fall within their windows. PFOA and
 670 PFOS must elute within 0.1 minutes of their associated EIS compounds.

671 **10.3** Initial Calibration672 **10.3.1** Initial Calibration Procedure

673 The most recent mass calibration must be used. Prior to initial calibration (ICAL), confirm the
 674 mass calibration used meets all manufacturer's criterion for confirmation of the mass calibration,
 675 mass resolution, and peak relative response. ICAL must be performed after a new mass
 676 calibration has been performed, whenever the laboratory has taken corrective actions that might
 677 affect the initial calibration criteria, and when the acceptance criterion of calibration verifications
 678 (ICV, CCV, and ISC) cannot be met. Following the analysis of an IB (Section 7.2.9), analyze
 679 each of the calibration standards that were prepared as described in Section 7.2.5. After the
 680 analysis of the highest level calibration standard, analyze another IB followed by an ICV. The
 681 same injection volume must be used for all standards, AFFF samples, blanks, and QC samples.

682 **10.3.2** Initial Calibration Calculations

683 PFOA and PFOS are quantified by isotope dilution quantitation, whereby the response of PFOA
 684 and PFOS are compared to the response of their associated EIS compounds. Calculate the
 685 response ratio (RR) for PFOA and PFOS in each calibration standard using the equation below.

$$686 \quad RR = \frac{Area_n M_l}{Area_l M_n}$$

687 where:

688 $Area_n$ = The measured area of the Q1 m/z for the native (unlabeled) PFAS

689 $Area_l$ = The measured area at the Q1 m/z for the corresponding isotopically labeled PFAS
 690 added to the sample before extraction

691 M_l = The mass of the isotopically labeled compound in the calibration standard

692 M_n = The mass of the native compound in the calibration standard

693

694 Calculate the response factor (RF_s) for each EIS in the calibration standard using the equation
 695 below.

$$696 \quad RF_s = \frac{Area_l M_{NIS}}{Area_{NIS} M_l}$$

697 where:

698 $Area_l$ = The measured area of the Q1 m/z for the isotopically labeled PFAS standard added
 699 to the sample before extraction

700 $Area_{NIS}$ = The measured area at the Q1 m/z for the isotopically labeled PFAS used as the non-
 701 extracted internal standard (NIS)

702 M_{NIS} = The mass of the isotopically labeled compound used as the non-extracted internal
 703 standard (NIS) in the calibration standard

704 M_i = The mass of the isotopically labeled PFAS standard added to the sample before
705 extraction

706 10.3.3 Instrument Linearity

707 10.3.3.1 To establish acceptable linearity, the calibration must meet one of the following
708 criteria:

709 Option 1: The calculated relative standard deviation (RSD) of the RR values of the
710 calibration standards for PFOA, PFOS and EIS and NIS compound each must be \leq
711 20%.

712 Option 2: The calculated relative standard error (RSE) of the calibration standards for
713 PFOA, PFOS and EIS and NIS compound each must be \leq 20%.

714 Option 3: Linear or non-linear regression calibration must have a coefficient of
715 determination (r^2) that is \geq 0.99 for PFOA and PFOS both.

716 10.3.3.2 In addition, to meeting the criteria stated in one of the options listed above, PFOA and
717 PFOS must recover within 70-130% of their true value for each calibration standard
718 included in the calibration curve. These criteria must be met before any AFFF samples,
719 QC samples, additional blanks, or calibration verifications are analyzed.

720 10.3.4 Initial Calibration Verification

721 10.3.4.1 IBs analyzed immediately before and after the initial calibration standards must not
722 contain PFOA or PFOS at a concentration \geq $\frac{1}{2}$ the LOQ in order the calibration to be
723 used. If the IB analyzed after the highest level calibration does not meet this criterion,
724 the instrument must be recalibrated using a series of calibration standards in which the
725 highest level standard is a lower concentration than the one previously used.

726 10.3.4.2 The ICV must meet the criteria stated in Section 9.2.1.2. If this criterion is not met,
727 corrective action must be taken and the ICV must be reanalyzed. If the ICV continues
728 to fail, initial calibration must be performed again.

729

730 11.0 Sample Preparation

731 11.1 Precautions

732 Fluoropolymer articles or *Kim-wipes* must not be used during sample preparation. Only HDPE or
733 polypropylene squeeze bottles and centrifuge tubes are to be used. Reagents and solvents for cleaning
734 syringes may be kept in glass containers.

735 11.2 Subsample Preparation

736 11.2.1 To prepare the sample, place a 60 mL PFAS-free HDPE sample bottle onto a four decimal place
737 balance and tare the weight. Place an aliquot of approximately 0.02 g of the AFFF sample into
738 the tared sample bottle and record the weight of the sample. This is the sample size.

739 11.2.2 Add 60 mL of PFAS-free reagent water to the bottle and swirl to mix.

740 11.2.3 To prepare the SD, repeat Section 11.2.1 and 11.2.2 using a second aliquot of the sample.

- 741 **11.2.4** Repeat the process in Sections 11.2.1 through 11.2.3 for each AFFF sample to be included in the
742 preparation batch (up to 20 AFFF samples).
- 743 **11.2.5** To prepare the MB, repeat Section 11.2.1 and 11.2.2 using a 0.02 g of PFAS-free reagent water as
744 the collected sample. Add 60 mL of PFAS-free reagent water to the bottle and swirl to mix.
- 745 **11.2.6** To prepare the LCS, LCSD, and LOQVER, for each, repeat the process in Section 11.2.1 using a
746 0.02 g of PFAS-free reagent water as the collected sample. Using a native spiking standard
747 solution, spike the LCS and LCSD at a concentration that is $> \text{LOQ}$ and \leq mid-level calibration
748 standard and the LOQVER at a concentration $\geq \text{LOQ}$ and ≤ 2 times the LOQ. Add 60 mL of
749 PFAS-free reagent water to the bottle and swirl to mix.
- 750 **11.2.7** Do not continue processing before at least 3 hours has lapsed to ensure the samples are
751 completely dissolved. Some samples may take additional time to dissolve. The dissolution time
752 for QC samples must be the same as the longest dissolution time required for an AFFF sample in
753 the batch. Inspect the samples to ensure that it has dissolved before proceeding. If it is not
754 completely dissolved, record this and notify the client before continuing. The 60 mL HDPE
755 sample bottle is now considered the AFFF sample container.
- 756 **11.2.8** After the required dissolution time has lapsed (Section 11.2.7), spike an aliquot of EIS solution
757 directly into each of the 60 mL AFFF sample containers. Swirl the samples to mix.
- 758 **11.2.9** Verify, using pH paper, that the pH of each sample is 6.5 ± 0.5 . Adjust the pH if necessary, with
759 50% formic acid, ammonium hydroxide, or with 5% formic acid and 3% aqueous ammonium
760 hydroxide.
- 761 **11.3** Sample Extraction and Cleanup
- 762 This section applies to all batch QC samples (MB, LCS, LCSD, LOQVER, and SDs) and AFFF samples.
763 Note that the volumes below are associated with validation of the specific cartridge used in validation.
764 150 mg WAX cartridges from other manufacturers may on different volumes. In addition, laboratories
765 may implement steps to reduce the final extract volume (e.g., by using nitrogen blowdown techniques) to
766 improve sensitivity. Addition of final extract volume reduction steps and changes to the volumes cited in
767 this section are permitted so long as all of the quality control criteria of this method are met.
- 768 **11.3.1** Collect the required number (one per sample) of Waters Oasis 150 mg WAX SPE cartridges.
769 Pack solvent rinsed silanized glass wool up to half the height of each cartridge barrel. Label each
770 cartridge with a sample identifier. Load the cartridges in the vacuum manifold and add a
771 reservoir and reservoir adaptor to each cartridge.
- 772 **11.3.2** To pre-condition the cartridges, without using the vacuum, wash the cartridges with 15 mL of 1%
773 methanolic ammonium hydroxide followed by 5 mL of 0.3 M formic acid. Do not allow the
774 cartridge to go dry. Discard the wash solvents.
- 775 **11.3.3** Pour each sample into its corresponding reservoir (do not use a pipette) while avoiding splashing
776 the sample. Empty the bottle as much as possible. Adjust the vacuum and pass the sample
777 through the cartridge at approximately 5 mL/min. Retain the emptied sample bottle and allow it
778 to air dry. These containers will be rinsed later (Section 11.3.5). Discard the sample being pulled
779 through the cartridge.
- 780 **11.3.4** Leaving the reservoir in place until after the sample has been eluted, rinse the walls of the
781 reservoir thoroughly with 5 mL reagent water twice (10 mL total) followed by 5 mL of 1:1 0.1M

782 formic acid:methanol and pass these rinses through the cartridge at a rate of approximately 5
783 mL/min using vacuum. Dry the cartridge by pulling air through for 15 seconds. Discard these
784 rinses.

785 **11.3.5** Label a 13 x 100 mm polypropylene collection tube (Section 6.9) for each sample and place them
786 in the manifold rack. Make sure extract delivery needles are positioned inside each collection
787 tube but are not touching the walls of the tubes. Rinse the inside of each of the 60 mL sample
788 bottles that were retained in Section 11.3.3 with a total of 5 mL of elution solvent (1% methanolic
789 ammonium hydroxide). Use a glass pipette to transfer the bottle rinse to its corresponding
790 reservoir, washing the walls of the reservoir. Use vacuum to pull the elution solvent through the
791 cartridges and into their corresponding collection tubes. After this step, the sample bottles no
792 longer need to be retained.

793
794 **11.3.6** Add 25 µL of acetic acid to each collection tube and vortex to mix. Using a 10 mg scoop, add 10
795 mg of carbon to a sample. Occasionally hand-shake the sample. Hand-shake for no longer than a
796 total of 5 minutes and immediately vortex for 30 seconds and centrifuge at 2800 rpm for 10
797 minutes.

798 **11.3.7** Label a fresh set of collection tubes. Use a 5 mL polypropylene syringe equipped with a syringe
799 filter (25 mm filter, 0.2 µm nylon membrane) to filter the entire extract into the prepared
800 collection tubes. Add NIS solution to each collection tube and vortex to mix. Transfer a portion
801 into a 1 mL polypropylene microvial for LC-MS/MS analysis. Cap the collection tube with the
802 remaining extract and store at 4°C.

803

804 **12.0 Instrumental Analysis**

805 The analysis of sample extracts for PFAS by LC-MS/MS is performed on a high performance or ultrahigh
806 performance liquid chromatograph (HPLC/UPLC) coupled to a triple quadrupole mass spectrometer,
807 running manufacturer's software. The mass spectrometer is run with unit mass resolution in the multiple
808 reaction monitoring (MRM) mode.

809 **12.1** Perform the mass calibration and mass calibration verification (Section 10.1), establish the
810 operating conditions (Section 10.2), and perform the initial calibration (Section 10.3). Samples can be
811 analyzed only after all associated criteria is met. The injection volume of samples (including batch QC
812 samples) must be the same as that of the standards and blanks. Standards and sample extracts must be
813 brought to room temperature and vortexed prior to aliquoting into an instrument vial in order to ensure
814 their homogeneity. Samples must be sequenced for analysis as follows

- 815 1. Instrument Blank (IB)
- 816 2. Instrument Sensitivity Check (ISC)
- 817 3. Continuing Calibration Verification (CCV)
- 818 4. PFOA Qualitative Identification Standard
- 819 5. Instrument Blank (IB)
- 820 6. Method Blank (MB)
- 821 7. Limit of Quantitation Verification Sample (LOQVER)
- 822 8. Laboratory Control Sample (LCS)
- 823 9. Laboratory Control Sample Duplicate (LCSD)
- 824 10. Samples (10 or fewer)

- 825 11. Sample Duplicates for samples analyzed in within bracket
826 12. Continuing Calibration Verification (CCV)
827 13. Instrument Blank
828 14. Samples (10 or fewer)
829 15. Sample Duplicates for samples analyzed in within bracket
830 16. Continuing Calibration Verification (CCV)
831 17. Instrument Blank (IB)

832
833 **12.2** All analytical batch QC standards (IB, ISC, and CCV) must meet the acceptance criteria stated in
834 Section 9.2.1 in order for sample results to be reported.

835

836 **13.0 Data Analysis, Calculations, and Reporting**

837 **13.1 Identification of Peaks**

838 AFFF samples may contain both branched and linear isomers of PFOA and PFOS, therefore all isomers
839 must be included in the quantitation of PFOA and PFOS. Retention times are established by the linear
840 isomer (Section 10.2.2.1), however, the retention time window for PFOS must encompass the retention
841 time of the branched isomers identified in the PFOS calibration standards while the retention time
842 window for PFOA must encompass the retention time of the branched isomers identified in the PFOA
843 qualitative standard analyzed after the initial calibration.

844 PFOA, PFOS, EIS compounds, and NIS compounds are considered identified in IBs, standards, AFFF
845 samples, and QC samples if the following criteria are met.

846 **13.1.1** The peak response must be at least three times the background noise level (S/N 3:1). This criteria
847 applies to both the quantification and confirmation ion peaks. If the S/N ratio is not met but the
848 background is low, then the analyte is to be considered a non-detect.

849 **13.1.2** AFFF samples may contain both branched and linear isomers of PFOA and PFOS, therefore all
850 isomers must be included in the quantification of PFOA and PFOS. The retention time window
851 for PFOS must encompass the retention time of the branched isomers identified in the PFOS
852 calibration standards while the retention time window for PFOA must encompass the retention
853 time of the branched isomers identified in the PFOA qualitative standard analyzed after the initial
854 calibration.

855 **13.1.3** The RTs of the PFOA, PFOS, EIS compound, and NIS compound peaks must fall within ± 0.4
856 minutes of the predicted RTs from the midpoint standard of the ICAL or initial daily CCV,
857 whichever was used to establish the RT window position for the analytical batch. In addition,
858 PFOA and PFOS peaks must elute within ± 0.1 minutes of the associated EIS. If this criteria is
859 not met, then the analyte is to be considered a non-detect.

860 **13.1.4** The ion ratio for each method analytes (PFOA and PFOS) in each sample (including preparatory
861 batch QC samples) must be determined as the ratio of the total (branched and linear isomers)
862 quantification ion response of the analyte to the total (branched and linear isomer) confirmation
863 ion response of the analyte. This ratio must fall within $\pm 50\%$ of the ion ratio observed in the mid-
864 level initial calibration standard. If the ion ratio failed to meet this criteria, reanalyze to confirm
865 the failure using a fresh aliquot of the extract. If the preparatory batch QC sample failure is

866 confirmed and ion ratios criteria is met for an associated sample, that sample may be reported.
 867 Confirming failures in preparatory batch QC samples and AFFF samples must report data with a
 868 “I” data qualifier and must be discussed in the case narrative.

869 13.2 Quantitative determination

870 Concentration of PFOA and PFOS are determined with respect to their corresponding EIS compounds.
 871 The EIS compound recoveries are determined with respect to their corresponding NIS compounds.

872 **13.2.1** Results for native compounds are recovery corrected by the method of quantification.

873 To calculate method analytes (PFOA and PFOS) concentrations:

$$874 \text{Concentration (ng/g)} = \frac{\text{Area}_n M_l}{\text{Area}_l(\overline{RR})} \times \frac{1}{W_s}$$

875 *where:*

876 Area_n = The measured area of the Q1 m/z for the native (unlabeled) PFAS

877 Area_l = The measured area at the Q1 m/z for the isotopically labeled PFAS (EIS).

878 M_l = The amount of the isotopically labeled compound added (ng)

879 \overline{RR} = Average response ratio used to quantify method analyte by the isotope dilution method

880 W_s = Sample weight (g)

881 **13.2.2** Extracted internal standard (EIS) recoveries are determined similarly against the non-extracted
 882 internal standard (NIS) and are used as general indicators of overall analytical quality.

883 To calculate EIS compounds ($^{13}\text{C}_8$ -PFOA and $^{13}\text{C}_8$ -PFOS) concentrations:

$$884 \text{Concentration (ng/g)} = \frac{\text{Area}_l M_{nis}}{\text{Area}_{nis} \overline{RF}_s} \times \frac{1}{W_s}$$

885 *where:*

886 Area_l = The measured area at the Q1 m/z for the isotopically labeled PFAS (EIS)

887 Area_{nis} = The measured area of the Q1 m/z for the non-extracted internal standard (NIS)

888 M_{nis} = The amount of the added non-extracted internal standard (NIS) compound (ng)

889 W_s = Sample weight (g)

890 \overline{RF}_s = Average response factor used to quantify the isotopically labeled compound

891 **13.2.3** If the response of a method analyte exceeds the quantitation range for any sample, extracts must
 892 be diluted to bring the exceeding analyte’s response within the calibration range, or the sample
 893 must be re-extracted for the exceeding analyte using a smaller aliquot of the AFFF sample.
 894 Instrument blanks must be analyzed following quantitation range exceedances per Section
 895 9.2.1.1.

896 13.3 Reporting

- 897 **13.3.1** Results must be reported in ppb. All associated QC data must be reported with the sample
898 results. Refer to the DoD/DOE QSM, Appendix A for general reporting requirements. Reports
899 must include all documentation needed to facilitate Stage 4 Validation per the *DoD Data*
900 *Validation Guidelines Module 3: Data Validation Procedure for Per- and Polyfluoroalkyl*
901 *Substances Analysis by QSM Table B-15* (Reference 17.3).
- 902 **13.3.2** Report concentrations \geq the DL to 3 significant figures. If the result falls between the DL and
903 LOQ, report the data with a “J” data qualifier, indicating it is an estimated value. If the method
904 analyte was not detected or the result falls below the DL, report the DL concentration with a “U”
905 data qualifier.
- 906 **13.3.3** The DL and LOQ values associated with each sample must take into account the volume of
907 sample prepared and any dilutions made. Report the result for each method analyte from the
908 lowest dilution that is within the quantification range that meets the acceptance criteria for EIS
909 and NIS recoveries. Samples with an associated LOQ that is \geq 25 ppb for PFOA or PFOS cannot
910 be used to demonstrate compliance to the requirements of MIL-PRF-24385.
- 911 **13.3.4** Results from tests performed with an analytical system that is not in control (i.e., that does not
912 meet acceptance criteria for any QC tests in this method) must be documented and reported (e.g.,
913 as a qualifier on results and discussed in the case narrative). Results with associated ion ratio,
914 EIS compound recovery, NIS compound recovery, and SD % RPD acceptance criteria
915 exceedances may not be used to demonstrate compliance to the requirements of MIL-PRF-24385.
916 Samples with an associated LOQ that is \geq 25 ppb for PFOA or PFOS cannot be used to
917 demonstrate compliance to the requirements of MIL-PRF-24385. All other failures must be
918 evaluated by NAVSEA to determine if the results are acceptable to demonstrate compliance.

919

920 **14.0 Method Performance**

921 Routine method performance must be monitored through certified reference materials (CRMs), when
922 commercially available. Currently, no such CRMs are available. When QC samples analyzed with
923 samples monitor ongoing method performance.

924 This method was validated using data generated by NAVSEA’s interlaboratory validation study.

925

926 **15.0 Pollution Prevention**

927 Application of this method must be compliant with all Federal, Provincial/State and Municipal
928 regulations governing waste management, including land disposal restrictions and sewage discharge
929 regulations. All standards are prepared in volumes consistent with volumes required by the method to
930 minimize the disposal of standards. The laboratory safety manual governs the safe storage, labeling and
931 disposal of laboratory wastes.

932

933 **16.0 Waste Management**

934 Application of this method must be compliant with all Federal, State, and local regulations governing
935 waste management, particularly the hazardous waste identification rules and land disposal restrictions,

936 and to protect the air, water, and land by minimizing and controlling all releases from fume hoods and
937 bench operations. Compliance is also required with any sewage discharge permits and regulations.

938

939 **17.0 References**

940 17.1 Department of Defense/Department of Energy Consolidated Quality Systems Manual for
941 Environmental Laboratories, <https://www.denix.osd.mil/edqw/>

942 17.2 Fire Extinguishing Agent, Aqueous Film-Forming Foam (AFFF) Liquid Concentrate, for Fresh
943 and Sea Water, MIL-PRF-24385, current version,

944

945 17.3 DoD Data Validation Guidelines Module 3: Data Validation Procedure for Per- and
946 Polyfluoroalkyl Substances Analysis by QSM Table B-15, <https://www.denix.osd.mil/edqw/>

DRAFT